

Quanterix

simoa

Direct Detection of Bacterial DNA and Viral RNA at Subfemtomolar Concentrations Using Single Molecule Arrays

Oak Ridge Conference
April 18, 2013

Agenda


- Introduction of single molecule arrays (Simoa)
- Assay process for Simoa DNA/RNA detection
 - Fragmentation of target DNA/RNA
 - Improved efficiency of capture and enzyme labeling of hybridized complexes on beads
 - Capture target DNA/RNA on paramagnetic beads
 - Hybridize biotinylated detection probes to captured targets
 - Label hybridized complexes with an enzyme
- Analytical sensitivity with purified DNA/RNA
- Detection of *S. aureus* in whole blood and environmental water

Counting Single Molecules Using Simoa

Rissin *et al.*, *Nat. Biotechnol.* **2010**, *28*, 595-599
Kan *et al.*, *Lab Chip* **2012**, *12*, 977-985

Quanterix Confidential | 3

Fragmentation of target genomic DNA/RNA




The samples were either fragmented at high concentrations and then diluted before testing in Simoa ("fragment-then-dilute", red squares), or diluted to the concentrations to be tested and fragmented ("dilute-then-fragment", blue circles).

Quanterix

Confidential | 7

Capture of fragmented DNA/RNA on beads - use of multiple capture probes




The fragmented DNA was captured by mixtures of multiple capture beads, each subpopulation of beads presenting a specific capture probe complementary to the different sites along the ~1,000 bp target gene. The number of beads was kept constant at 500,000 beads per 100 mL sample.

Quanterix

Confidential | 8

Labeling of hybridized complexes - use of multiple detection probes




- Single biotin vs. multiple biotins on detection probes
 square: (B)iotin-aaagaagaggtgttagttatgac
 triangle: (B)iotin-aaagaagaggt(B)gttagtt(B)atgac
 circle: (B)iotin-aaagaagaggtgttagttatgac; (B)-taagtgtgcatatgtatgac;
 (B)-aggatggctatcagtaattt

Quanterix


Confidential | 9

Labeling of hybridized complexes
 - use of optimized enzyme concentration




The highest S/B ratio (16.3) was observed at [SβG] = 200 pM.

Dose-response curve
 - detection of purified genomic DNA from *S. aureus*



•LOD (3xSD+ background): **0.04fM**, equivalent to **1200 DNA molecules/50 μl** sample volume, for this run.

Dose-response curve
 - detection of purified genomic RNA from Sendai virus



•LOD (3xSD+ background): **0.06fM**, equivalent to **1800 RNA molecules/50 μl** sample volume, for this run.
