



**Next Generation Protein Detection**

Oak Ridge Conference  
Introducing Longitudinal Assay Screening (LAS)  
April 19, 2013


---

---

---

---

---

---

---

---

**Introduction to Inanovate**

<b>Location</b>	Research Triangle, NC; Boston, MA
<b>Technology</b>	<b>Longitudinal Assay Screening</b>
<b>Associated Product</b>	Bio-ID 400
<b>Application Focus</b>	Detection, measurement and monitoring of multiple proteins
<b>Under Development</b>	<ul style="list-style-type: none"> <li>Automated near-care platform (Bio-ID Dx)</li> <li>Diagnostic biomarker panels for prostate cancer, ovarian cancer and sepsis.</li> </ul>
<b>Selected Partnerships &amp; Collaborations (public)</b>	MD Anderson, Harvard Medical School, Thermo Fisher, Dana Faber Cancer Center, Manchester Hospital (NHS Trust)

Proprietary and Confidential 2

---

---

---

---

---

---

---

---

**Protein Multiplexing: The Problems**

Accurately detecting and measuring single, let alone **multiple proteins** is a significant scientific and technical achievement.

The development of Longitudinal Assay Screening (LAS) was driven by three limitations with existing technologies

Limited Detection Range

Limited Biological Relevance

Limited Data Accuracy

Proprietary and Confidential 3

---

---

---

---

---

---

---

---

## Protein Multiplexing: The Problems

Image borrowed from Journal of Proteome Research • 2011, 10, 5-16 • Sibihi Sathirava et al

Development of Plasma Protein Biomarkers




Figure 1. Dynamic plasma protein concentration range and the dynamic range of plasma protein categories are shown as reported by Anderson et al.<sup>16</sup> Red dots indicate proteins identified by the HUPO plasma proteome project (PP1)<sup>16</sup> and red line represent commonly used biomarkers in the clinic. Suitable minimal range of biomarker screening in plasma is shown with dotted lines. Adapted from Sathirava et al.<sup>16</sup>

Limited Detection Range

Auto-immune diseases, e.g. rheumatoid arthritis.  
Inflammatory markers: Impact cardiovascular disease to cancer.

Limited Biological Relevance

Proprietary and Confidential

4

---

---

---

---

---

---

---


---

---

---

## Protein Multiplexing: The Problems

Which signal is 'real' and which is non-specific background?



Limited Data Accuracy

Increases time and cost of biomarker discovery and validation programs, and introduces 'false positive errors' into clinical applications.

Proprietary and Confidential

5

---

---

---

---

---

---

---

---

---

---

## The Inanovate Solution: Longitudinal Assay Screening (LAS)

Inanovate has overcome the problems of protein multiplexing through the development of a new category of protein screening technology: Longitudinal Assay Screening (LAS)

The Bio-ID



A new solution to protein detection

Inanovate has recently completed testing and demonstration of the world's first protein detection platform to integrate LAS technology: The Bio-ID 400.

Proprietary and Confidential

6

---

---

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

### The Bio-ID: First platform to integrate LAS technology

**The Bio-ID Detector**

The hardware that detects, processes and analyzes the assays on the Bio-ID's disposable 'Chips'

**The Bio-ID Disposable 'Chip'**

A fluidic cartridge that houses the assays and facilitates compatibility with the Detector.

**Real-time Biomarker Measurement**

Iterative flow of sample and detection antibody across the surface of the protein array

Real-time detection, kinetic signatures, increased confidence, walkaway automation!

VIDEO

Proprietary and Confidential 8

---

---

---

---

---

---

---

---

### Analyzing LAS data

The Bio-ID captures time course fluorescent data

Rates are calculated and plotted vs. concentration

The rate based binding curve is used to quantify unknowns

IL6	
Vmax	1255.10
n	0.974
Km	2568.39
HOD	53.233
LDD	0.85
LQD	1.70

Proprietary and Confidential 9

---

---

---

---

---

---

---

---

## The Benefits

- Rate based analysis enables accurate quantitation across over a 7log range in a single multiplex. Eliminates need for serial dilutions, making multiplexing faster, cheaper and helping preserve precious samples.
- Large detection range, LAS allows users to run virtually any assay of interest in one test, enabling the development of truly biologically relevant multiplexes.
- Analysis of real-time kinetic data improves identification and discrimination of background and non-specific signals, delivering more accurate quantitation at low analyte concentrations.

Limited Detection Range

Limited Biological Relevance

Limited Data Accuracy

Proprietary and Confidential 10

---

---

---

---

---

---

---

---

---

---

## Demonstration data: Improving detection range & biological relevance

To illustrate the Bio-ID's capacity to quantitate very different protein concentrations in a single run, known low abundance proteins (IL-6 and IL-1b) were analyzed in the same 3-plex assay as a known high abundance protein (CRP)

Assay	LDD	HQD	Inter-run %CV (cross instrument)
IL-6	0.85	53,000	3.0
IL-1b	1.26	21,000	6.0
CRP	21.89	1,361,000	5.3

The Bio-ID 400 is an early stage technology. Significant further improvement is likely as the system and assay are further optimized through the coming months.

**Large detection range = resource savings + biologically relevant assays**

Proprietary and Confidential 11

---

---

---

---

---

---

---

---

---

---

## Demonstration data: Improving detection range & biological relevance

Despite CRP being present at up to 1,000,000X the concentration of IL-6 and IL-1b, there is no significant cross-talk from CRP to the IL-6 and IL-1b assays.

CRP Crosstalk

IL6 Crosstalk

IL1b Crosstalk

Proprietary and Confidential 12

---

---

---

---

---

---

---

---

---


---



**Summary & Conclusions:  
Demo 3-plex data from the Bio-ID 400**

Assay	LDD	HQD	Inter-run %CV	% Recovery
IL-6	0.85	53,000	3.0	96.7
IL-1b	1.26	21,000	6.0	97.6
CRP	21.89	1,361,000	5.3	94.4

Large detection range = resource savings + biologically relevant assays



Kinetic analysis = high quality data, non-targeted discrimination & improved accuracy

Proprietary and Confidential

16

---

---

---

---

---

---

---

---

---

---

**Thank you for your attention**



**Contact Information:**

David Sloan  
Director of Assay Development  
[dsloan@inanutate.com](mailto:dsloan@inanutate.com)  
919-249-8637

David Ure  
President and CEO  
[dure@inanutate.com](mailto:dure@inanutate.com)  
617-610-1712

---

---

---

---

---

---

---

---

---

---