



# Next Generation Protein Detection

**Oak Ridge Conference**  
**Introducing Longitudinal Assay Screening (LAS)**

*April 19, 2013*



**inanovate**  
Opening doors to a smaller world



# Introduction to Inanovate



<b>Location</b>	Research Triangle, NC; Boston, MA
<b>Technology</b>	<b>Longitudinal Assay Screening</b>
<b>Associated Product</b>	Bio-ID 400
<b>Application Focus</b>	Detection, measurement and monitoring of multiple proteins
<b>Under Development</b>	<ul style="list-style-type: none"><li>• Automated near-care platform (Bio-ID Dx)</li><li>• Diagnostic biomarker panels for prostate cancer, ovarian cancer and sepsis.</li></ul>
<b>Selected Partnerships &amp; Collaborations (public)</b>	MD Anderson, Harvard Medical School, Thermo Fisher, Dana Faber Cancer Center, Manchester Hospital (NHS Trust)



# Protein Multiplexing: The Problems

Accurately detecting and measuring single, let alone **multiple proteins** is a significant scientific and technical achievement.

The development of Longitudinal Assay Screening (LAS) was driven by three limitations with existing technologies

**Limited Detection  
Range**

**Limited Biological  
Relevance**


**Limited Data  
Accuracy**

# Protein Multiplexing: The Problems

Image borrowed from Journal of Proteome Research - 2011, 10, 5-16 - Silvia Surinova et al

Development of Plasma Protein Biomarkers

reviews




**Figure 1.** Dynamic plasma protein concentration range and the three main plasma protein categories are shown as reported by Anderson et al.<sup>18</sup> Red dots indicate proteins identified by the HUPPO plasma proteome project (PPP)<sup>109</sup> and yellow dots represent currently used biomarkers in the clinic.<sup>6</sup> Suitable minimal range of detection for biomarker targeting in plasma is shown with dotted lines. Adapted from Schiess et al.<sup>4</sup>

**Limited Detection Range**

**Auto-immune diseases, e.g. rheumatoid arthritis.**

**Inflammatory markers: Impact cardiovascular disease to cancer.**


**Limited Biological Relevance**



# Protein Multiplexing: The Problems



Which signal is 'real' and which is non-specific background?



Limited Data Accuracy


Increases time and cost of biomarker discovery and validation programs, and introduces 'false positive errors' into clinical applications.

# The Inanovate Solution: Longitudinal Assay Screening (LAS)

Inanovate has overcome the problems of protein multiplexing through the development of a new category of protein screening technology:

**Longitudinal Assay Screening (LAS)**


## The Bio-ID



**A new solution  
to protein  
detection**

Inanovate has recently completed testing and demonstration of the world's first protein detection platform to integrate LAS technology: **The Bio-ID 400**.

# What is LAS?



# The Bio-ID: First platform to integrate LAS technology

## The Bio-ID Detector




The hardware that detects, processes and analyzes the assays on the Bio-ID's disposable 'Chips'

## The Bio-ID Disposable 'Chip'



A fluidic cartridge that houses the assays and facilitates compatibility with the Detector.

## Real-time Biomarker Measurement



Iterative flow of sample and detection antibody across the surface of the protein array

Real-time detection, kinetic signatures, increased confidence, walkaway automation!


**VIDEO**

Proprietary and Confidential




# Analyzing LAS data

The Bio-ID captures time course fluorescent data



Rates are calculated and plotted vs. concentration



The rate based binding curve is used to quantify unknowns

**IL6**

<b>Vmax</b>	1255.10
<b>n</b>	0.974
<b>Km</b>	2588.39

<b>HQD</b>	53,233
<b>LDD</b>	0.85
<b>LQD</b>	1.70



# The Benefits

- Rate based analysis enables accurate quantitation across over a 7log range in a single multiplex. Eliminates need for serial dilutions, making multiplexing faster, cheaper and helping preserve precious samples.
- Large detection range, LAS allows users to run virtually any assay of interest in one test, enabling the development of truly biologically relevant multiplexes.
- Analysis of real-time kinetic data improves identification and discrimination of background and non-specific signals, delivering more accurate quantitation at low analyte concentrations.

**Limited Detection Range**

**Limited Biological Relevance**

**Limited Data Accuracy**



# Demonstration data: Improving detection range & biological relevance

To illustrate the Bio-ID's capacity to quantitate very different protein concentrations in a single run, known low abundance proteins (IL-6 and IL-1b) were analyzed in the same 3-plex assay as a known high abundance protein (CRP)

Assay	LDD	HQD	Inter-run %CV (cross instrument)
IL-6	0.85	53,000	3.0
IL-1b	1.26	21,000	6.0
CRP	21.89	1,361,000	5.3




The Bio-ID 400 is an early stage technology. Significant further improvement is likely as the system and assay are further optimized through the coming months.

**Large detection range = resource savings + biologically relevant assays**



# Demonstration data: Improving detection range & biological relevance

Despite CRP being present at up to 1,000,000X the concentration of IL-6 and IL-1b, there is **no significant cross-talk from CRP to the IL-6 and IL-1b assays.**






# Demonstration data:

## Improving detection range & biological relevance


Assay	Level	% Recovery	
		Simulated Matrix	Human Serum
IL-6	High	93.2	108.8
	Mid	75.4	72.2
	Low	96.9	117.3
IL-1b	High	97.6	102.6
	Mid	79.9	99.2
	Low	94.7	106.6
CRP	High	129.6	80.2
	Mid	94.4	90.2
	Low	88.9	73.4

The Bio-ID 400 is an early stage technology. Significant further improvement is likely as the system and assay are further optimized through the coming months.

# Demonstration data: Improving data accuracy



Orange data from known 30K pg/mL IL1-b standard. Blue from printed IL1-b detector Ab that reacts with streptavidin (simulated non-targeted signal).



## Assay Specific Signatures:

- Specific vs. non specific binding
- Curve shape
- Vmax
- Assay time-course



# Summary & Conclusions: The Benefits

- Rate based analysis enables accurate quantitation across over a 7log range in a single multiplex. Eliminates need for serial dilutions, making multiplexing faster, cheaper and helping preserve precious samples.
- Large detection range, LAS allows users to run virtually any assay of interest in one test, enabling the development of truly biologically relevant multiplexes.
- Analysis of real-time kinetic data improves identification and discrimination of background and non-specific signals, delivering more accurate quantitation at low analyte concentrations.

**Limited Detection  
Range**

**Limited Biological  
Relevance**



**Limited Data  
Accuracy**

# Summary & Conclusions:

## Demo 3-plex data from the Bio-ID 400

Assay	LDD	HQD	Inter-run %CV	% Recovery
IL-6	0.85	53,000	3.0	96.7
IL-1b	1.26	21,000	6.0	97.6
CRP	21.89	1,361,000	5.3	94.4

**Large detection range = resource savings + biologically relevant assays**



**Kinetic analysis = high quality data, non-targeted discrimination & improved accuracy**





Thank you for your attention



**Contact Information:**

David Sloan  
Director of Assay Development  
[dsloan@inanovate.com](mailto:dsloan@inanovate.com)  
919-249-8637

David Ure  
President and CEO  
[dure@inanovate.com](mailto:dure@inanovate.com)  
617-610-1712