

T4 CDC laboratory protocol for detecting carbapenem-resistant or carbapenemase-producing *Klebsiella* spp. and *E. coli* from rectal swabs

<p>STEP 1 Day One</p>	<p>Aseptically, place one 10-μg ertapenem or meropenem disc in 5 ml trypticase soy broth (TSB).</p> <p>Immediately inoculate the broth with the rectal swab.</p> <p>Incubate overnight at $35 \pm 2^\circ\text{C}$, ambient air.</p>
<p>STEP 2 Day Two</p>	<p>Vortex incubated broth culture and subculture 100 μl onto a MacConkey agar plate and streak for isolation.</p> <p>Incubate overnight at $35 \pm 2^\circ\text{C}$, ambient air.</p>
<p>STEP 3 Day Three</p>	<p>Examine the MacConkey agar for lactose-fermenting (pink-red) colonies. More than one colony morphology may represent different species of Enterobacteriaceae.</p> <p>It may be necessary to subculture representative colonies of each morphology type to non-selective media for isolation and/or for susceptibility testing.</p> <p>Screen representative isolated colonies using a phenotypic test for carbapenemase production, such as the modified Hodge test (MHT) or test for carbapenem susceptibility using a standardized method and follow the CLSI guidelines for identification of carbapenemase-producing Enterobacteriaceae.</p>
<p>STEP 4 Day Four</p>	<p>For CRE and/or MHT-positive isolates, perform species-level identification.</p>

Source: 2012 CRE Toolkit—Guidance for Control of Carbapenem-resistant Enterobacteriaceae (CRE)