

# PUTTING IT ALL TOGETHER: MOVING DATA AND REPORTING RESULTS

Judy Stone, Ph.D., MT (ASCP), DABCC  
Center for Advanced Laboratory Medicine  
University of California, San Diego  
jastone@ucsd.edu

---

---

---

---

---

---

---

## Financial disclosure

- Nothing to disclose

---

---

---

---

---

---

---

## Learning objectives

- Define the IT related terms commonly used in a clinical instrument interfacing project
- Explain the advantages of an instrument interface using ASTM/HL7 protocol as compared to a flat file export
- List 4 quality parameters used in LC-MSMS data review that can be automated

---

---

---

---

---

---

---

## Acknowledgements

- Tony DaSilva, Jonathan Agustin, Julia Drees at Kaiser Regional Laboratory
- Shamim Haider, Dave Lavorato, Chris Borton, Hua-fen Liu, John Gibbons at AB Sciex
- Mark George, Don Mason at Waters Corporation
- Don Jarvis at Data Innovations
- David Pope, Melodie Bayongan, John Gaede, Rob Fitzgerald at UC San Diego
- Mike Vogelsberg at SCC Soft Computer

---

---

---

---

---

---

---

## Outline

- How does interfacing your LC-MSMS differ from interfacing an automated chemistry analyzer? (*homework before you start the project*)
- Step by step – project planning to interface your LC-MSMS
- Automating LC-MSMS data review

---

---

---

---

---

---

---

## Acronyms

- Laboratory Information System (LIS)
- Laboratory Information Management System (LIMS)
- **LIS ≠ LIMS**
- Middleware (MW)

---

---

---

---

---

---

---

## The problem

- “Interface” is defined differently by

1. Clinical chemists, any other laboratorian
2. Laboratory LIS coordinators and LIS vendors
3. LC-MSMS vendors



---

---

---

---

---

---

---

---

## Common LIS Coordinators (Lab) and LIS vendors View

- Interface = connection of a medical device to the LIS for transfer of patient data using standardized protocols that are compliant with health information transfer regulations
- **RS232** connections - instrument to hospital LAN
- **ASTM** (or **HL7**) – standardized data exchange protocol (format & mapping)
- The instrument vendor and the LIS vendor will speak similar jargon, have similar expectations of the task

**Jargon Alert:** LAN = Local Area Network (Intranet)

---

---

---

---

---

---

---

---

**Standardized**

**Limited options**

**Regulated**

---

---

---

---

---

---

---

---

## Common LC-MSMS Vendors View

- Interface = export of LC-MSMS data in a **flat file** (generic format) to a shared network directory
- **TCP-IP** connections – instrument to hospital LAN
- .txt or .csv file format with (perhaps) a template customized by the LC-MSMS vendor to control what data fields are exported
- Purchase of an instrument does not automatically include responsibility of the LC-MSMS vendor to customize connectivity and data transfer to the LIS

---

---

---

---

---

---

---

**Jargon Alert:** flat file = A flat file is a set of data in a plain text file (.txt) . Each string in the text file holds one record, with fields separated by delimiters, such as commas or tabs.

**Jargon Alert:** TCP-IP = Transmission Control Protocol/Internet Protocol (*TCP/IP*) is the suite of communications protocols used to connect hosts (computers) on the Internet

---

---

---

---

---

---

---

## Generic

(.txt or .csv file format)



## Customizable



## Flexible



---

---


---

---

---

---

---



---

---


---

---

---

---

---



Common laboratorian view

• **“Not in my job description!”**

---

---


---

---

---

---

---



## Explaining the differences - 1

### LC-MSMS

- TCP-IP connections – instrument to hospital LAN

### Automated Chemistry Analyzer

- RS232 connections - instrument to hospital LAN

---

---

---

---

---


---

---

## Decoding the jargon-Connectivity Chemistry instrument connectivity

- RS232 - serial or COM port

(In Control Panel)



- Older
- Simpler
- Limited in distance from instrument to PC
- Must convert to TCP/IP for LAN connection

---

---

---

---

---

---

---

---

## Decoding the jargon-Connectivity

### LC-MSMS connectivity

- TCP/IP - network card (NIC) & ethernet cable



- Newer
- Direct connection to LAN, no conversion needed

---

---

---

---


---

---

---

---

## Decoding the jargon-Connectivity




---

---

---

---


---

---

---

---

## Decoding the jargon-Connectivity




---

---

---

---

---

---

---

---

## Connectivity barriers

- Your Network (LAN) Administrator requires automatic Windows Updates on the PC that controls the MS
- The MS vendor states “No Windows Updates”
- Your Network Administrator requires real time Antivirus scans
- The MS vendor states “Can’t install or run Antivirus”

---

---

---

---

---

---

---

---

## Work arounds - Antivirus

- Discuss with the MS vendor in advance (include in the RFP) & ask Network Admin. to be flexible about which Antivirus vendor
- Schedule Antivirus updates & PC scanning for times when the instrument is not acquiring
- You can exclude some directories of the MS vendor software from scanning – but that increases risk
- If possible – restrict internet sites on the MS PC

---

---

---

---

---

---

---

---

## Work arounds – OS Updates

- If a Windows update causes a problem – use Windows Restore Point function to return the PC to it's status quo before updates
- PC configuration - new restore point created before any Windows updates are installed
- Install only "Important", i.e. security related Windows updates
- Configure the PC to download & install, or download with manual install, some time after the 2<sup>nd</sup> Tues of each month whenever the instrument is not acquiring (e.g., Sunday 5am)

---

---

---

---

---

---

---

---

## Ease of access versus Security

- Security is a shared responsibility:  
IT, Laboratory, MS vendor
- Include IT in the planning for MS installation
  1. Get agreement on OS Updates & Antivirus
  2. Schedule activation of a Loc Plate near your instrument
  3. Obtain a **static IP address** from Network Admin. or schedule IT to configure for **dynamic IP address (DHCP)**
  4. **Vendor remote access to the MS should be controlled by the laboratory, NOT by the vendor (require "challenged access")**



**Jargon Alert: IP Address = PC (network card) unique identifier on the network 10.30.2.181**




---

---

---

---


---

---

---

---

## Instrument Network Architecture




---

---

---

---

---

---

---

---



## Explaining the differences - 1

### LC-MSMS

- **TCP-IP** connections – instrument to hospital LAN

### Automated Chemistry Analyzer

- **RS232** connections - instrument to hospital LAN



---

---

---

---

---

---

---

## Explaining the differences - 2

### LC-MSMS

- .txt or .csv file format with (perhaps) a template customized by the LC-MSMS vendor to control what data fields are exported

### Automated Chemistry Analyzer

- **ASTM** (or **HL7**) – standardized data exchange protocol (format & mapping)

---

---

---

---

---

---

---

## Why do standards & format matter?

- **Translation** – you call it MRM ratio and I call it 1° Target ratio and she calls it qualifier ratio
- **Location** = the address. Whatever you want to call it – how do I find it in the file?

The same logic for locating houses on a map:

- **Field Name/Tag/HL7 Segment** = apt #, street address, city, state, zip code
- **Field value** = John Doe  
(who lives there)



---

---

---

---


---

---

---

## Decoding the jargon – Data Exchange Standards

- **ASTM and HL7 are medical data exchange standards for all types of medical information (ASTM came before HL7)**



D Nelson, Development and Current Use of ASTM/HL7 Messages for Health Care Communications in the USA, *Stud Health Technol Inform*, 1993;6:138-43.

## Good news... but you have to ask:

- The LC-MS vendor may have an interface driver for your LIS/MW (**plug & play**) – ask!
- The LC-MS vendor may have a generic HL7 interface available – ask!
- The LIS vendor may have a robust utility for importing flat files (.txt, .csv) that have been exported from an instrument – ask!
- But be prepared for – flat file export, \$\$ and time for LIS vendor\* to build the interface

\*alternative: In-house programming expertise

Remember –  
your LIS  
coordinators  
and LIS vendor  
“speak” HL7!

## Why a flat file?

```
%PDF-1.3
%
1 0 obj
<<
/Linearized 1
/O 55
/H [ 2403 456 ]
/L 300652
/E 263086
/N 4
/T 299474
>>
endobj
```

```

Dlãj±á>þý
‘9;.<=>?@
ABCDEFGH-
!"#$$%&'()*+,-
./0123456789:;
<=>?@A B C D E F G H I J K L M N O P Q R S T U V W X Y Z
```

translation you call it mrm ratio and i call it 1\* target ratio and she calls it qualifier ratio location the address whatever you want to call it how do I find it in the text string the same logic for locating houses on a map field name tag hl7 segment apt street address city state zip code field value john doe who lives there

Convert formatted & structured to → generic

## .txt file output of an LC-MSMS result table

```
<?xml version="1.0"?>
<QANDASET description="" version="1"><XMLFILE filename="c:\users\jud\2c
modifiedtime="10:34:51"/><DATASET filename="c:\MasslynxProject\2013\0213.P4
creationdate="04 Feb 2013" creationtime="11:04:18"/><GROUPDATA count="1"><
MasslynxProject\2013\0213.P4\Methob\vitamin0_1f_2_no_smooth.mdb" modified
creationtime="10:28:10"/><SAMPLELISTDATA filename="c:\MasslynxProject\2013\
count="79" SAMPLE ID="1" groupid="1" name="02012013.1" createdate="01-Feb-
dilutionfactor="0.0000000000" extractvolume="0.0000000000" initamount="0.0000
stdconc="0.0000000000" stockdilutionfactor="0.0000000000" subjecttext="" subje
inletmethodname="c:\MasslynxProject\2013\0213.P4\ACQU08\vitamin04" msmethc
peakname="0.0000000000" peakid="0.0000000000" hplcmethodname="" tunemethodname="
instrument="Xevo-TQSGAA446" lab="" conditions="" submitter="" task="" user
type="" cas="" stdconc="0.0000000000"><PEAK foundscan="0" foundrt="0.000000
area="0.0000000000" height="0.0000000000" response="0.0000000000" pkflags="
concal="Nan" modifieddate="" modifiedtime="" modifiedtext="" modifieduser
starttime="0.0000000000" endtime="0.0000000000" absresponse="0.0000000000" rtrr
quarationwin="0.0000000000" ionratio="0.0000000000" ionratiopred="0.0000000
detectionthreshold="0.0000000000" detectionflag="0" quantthreshold="0.000000
chromtrace="395.400 &gt; 377.37" peaks="0" pkwidth="0.0000000000" pksigma="
heightdivarea="0.0000000000" baselinewidth="0.0000000000" peakquality="N/A
recoveryflag="0" matrixblankflag="0" solventblankflag="0" derivflag="0" deriv1
qcrtimeflag="0" qcpeakshapeflag="0" signoise="0" signoiseflag="0" cdflag="C
peakcalwidth="0.0000000000" peakasymmetryvalue="0.0000000000" peakrecovery
chromiseght="13950.0000000000" concdevpc="0.0000000000" lowerbound="0.
lowerbound="0.0000000000" medlumbound="0.0000000000" medlumbound2="0.0000
upperbound="0.0000000000" upperbound="0.0000000000" upperbound3="0.000000
toxconc1="0.0000000000" toxconc2="0.0000000000" toxconc3="0.0000000000" tox
toxfactor="0.0000000000" toxfactor4="0.0000000000" toxlod="0.0000000000"
toxlost="0.0000000000" toxlost2="0.0000000000" toxloq3="0.0000000000" toxloq
p1csreverse="0" ifit="N/A" ifitnorm="N/A" ifitconfidence="N/A" foundmass="N
ifitconfFlag="0" moerrorflag="0" gpcerrorflag="0"><CONFIRMATIONPEAK ic
predrt="0.0000000000" predrrt="0.9894999862" area="0.0000000000" height="0.
empc="0.0000000000" usanalconc="0.0000000000" concal="Nan" modifieddate="
starttime="0.0000000000" endtime="0.0000000000" starttime="0.0000000000" endtime="0.0000000000"></PEAK></GROUPDATA></QANDASET>
```

## More on file formats

• .txt or .csv flat file, most generic

• .xml “markup” language??

### XML:

- Data + Format rather than Data only
- Widely used
- Ultimately easier if people and software at both vendors can process XML

- XML (Extensible Markup Language) is a flexible way to create common information formats and share both the format and the data on the World Wide Web
- a set of rules for encoding documents in a format that is both human-readable and machine-readable
- is used to define documents with a standard format that can be read by any XML-compatible application. Basically, XML allows you to create a database of information without having an actual database.

## Flat file conversion to LIS format

1. Contract with the LIS vendor to do it OR
2. Use in-house programming resources to do it

a. Blank GE, Virji MA. *Development and implementation of an electronic interface for complex clinical laboratory instruments without a vendor-provided data transfer interface.* **J Pathol Inform** 2011;2:14

b. French D, Terrazas E. *The successful implementation of a licensed data management interface between a Sunquest® laboratory information system and an AB SCIEX TM mass spectrometer.* **J Pathol Inform** 2013;4:1

Convert generic to → formatted & structured

## Explaining the differences - 2

### LC-MSMS

- .txt or .csv file format with (perhaps) a template customized by the LC-MSMS vendor to control what data fields are exported

### Automated Chemistry Analyzer

- **ASTM** (or **HL7**) – standardized data exchange protocol (format & mapping)



---

---


---

---

---

---

---



---

---

---

---

---

---

---

### Interfacing Workflow A (results upload only)

1. Identify key technical contact people in your LIS and IT groups and at the LC-MS vendor
2. Gather and organize key information, set and communicate goals
  - Your IT group
  - LC-MS vendor
  - Your LIS group & LIS vendor
3. Kick off conference call with all parties (write up and distribute minutes)

---

---

---

---

---

---

---

## Interfacing Workflow B

4. Connect the MS to the network
5. Setup and test connectivity to the LIS/MW
6. Configure/Customize the data that will be exported from the MS  
(select data & format)
7. Mapping (translation tables)
8. Test Transmission in TEST LIS/MW system  
(send & look at communication traces)

---

---

---

---

---

---

---

## Interfacing Workflow C

9. Troubleshoot & repeat transmission
10. Decision – use LC-MSMS vendor, MW or LIS rules?
11. Test rules, fix rules, retest
12. Test batches, scale up in TEST system
13. Setup in LIVE, repeat all testing
14. Go live
15. Post-live monitoring



---

---

---

---

---

---

---

## Workflow detail

1. Gather and organize key information
  - a. IT
    - – what steps are needed to connect to your LAN (OS updates, Antivirus, type of IP address, restricting websites)?
    - Options for instrument network architecture (acquisition to a server)?
  - b. LC-MS vendor
    - – what interface options are there (LIS/MW driver, HL7 interface, simple flat file export)?
    - OS updates and antivirus compatibility?
    - Options for data delocalization (stand alone data analysis PCs or acquisition to a server)?
  - c. LIS – LIS vendor experience with LC-MS interfacing?

---

---

---

---

---

---

---

## Workflow detail

### 2. Connecting the MS to the network

- Two NIC cards in the instrument PC (one to the LC-MSMS and one to the LAN)
- Live loc plate & a network cable & an IP address
- TCP/IP connection – not RS232 connection
- Operating System (OS) i.e. Windows updates
- Antivirus installation & scanning



---

---

---

---

---

---

---

---

## Workflow Detail

### 3. Setup and test connectivity to the LIS/MW

Address for where to send the instrument data?

- Simple file export** – shared directory on the network
  - Copy/Paste a file to the shared directory  
(file permissions for writing a file to the shared directory?)
  - Can the LIS/MW capture and read the file?
- LIS/MW direct connection** – how to connect?
  - obtain the IP address & port # of LIS/MW
  - Enter them in to the instrument interface configuration screen  
(port # is virtual, think TV channel, not a physical plug-in port)
  - “Ping” the LIS/MW - does it “hear” you?

---

---

---

---

---

---

---

---

## Workflow Detail

### 4. How to find the IP address

- Click on  type in “cmd”

```
Ethernet adapter Instrument LAN:
Connection-specific DNS Suffix . :
Description . . . . . : Broadcom NetXtreme Gigabit Ethernet #2
Physical Address. . . . . : 00-52-8F-3D-83-81
Dhcp Enabled. . . . . : No
Autotest operation Enabled . . . : Yes
IP address. . . . . : 10.10.10.1 (Preferred)
Subnet Mask . . . . . : 255.255.255.0
Default Gateway . . . . . : 192.168.8.1 (Preferred)
Subnet Mask . . . . . : 255.255.255.0
Default Gateway . . . . . :
NetBIOS over Tcpip. . . . . : Disabled

Ethernet adapter Local Area Connection:
Connection-specific DNS Suffix . :
Description . . . . . : Realtek PCIe GbE Family Controller
Physical Address. . . . . : 00-52-8F-3D-83-81
Dhcp Enabled. . . . . : No
Autotest operation Enabled . . . : Yes
IP address. . . . . : 10.10.10.1 (Preferred)
Subnet Mask . . . . . : 255.255.255.0
Default Gateway . . . . . : 192.168.8.1 (Preferred)
Subnet Mask . . . . . : 255.255.255.0
```

### 5. How to ping the LIS/MW

- Go to CMD prompt ,
- type “ping ip address”  
(ping 10.10.10.1)

```
C:\Users\user>ping 10.10.10.1
Pinging 10.10.10.1 with 32 bytes of data:
Reply from 10.10.10.1: bytes=32 time=1ms TTL=64
Reply from 10.10.10.1: bytes=32 time=1ms TTL=64
Reply from 10.10.10.1: bytes=32 time=1ms TTL=64
Reply from 10.10.10.1: bytes=32 time=1ms TTL=64

Ping statistics for 10.10.10.1:
    Packets: Sent = 4, Received = 4, Loss = 0%
    Approximate round trip times in milliseconds:
        Minimum = 0ms, Maximum = 0ms, Average = 0ms
```

---

---

---

---

---

---

---

---

## Workflow detail

6. Configure/Customize export content & format from the MSMS



- What fields from your MS data to send (may be >100 fields that can be sent)?
- What are the field names?
- Will you send just the basics or more?



**Basics:**  
**Accession #**  
**Test name**  
**Test result**  
**Instrument ID**




---

---

---

---

---

---

---

---

## What is “more”?

Yes or No (Fail or Pass) responses for LC & MSMS quality parameters

- MRM qualifier ratio
- Expected retention time
- Relative retention time
- Internal standard peak area relative to calibrators
- Peak width
- Peak asymmetry
- Calibration characteristics (% deviation,  $R^2$ )
- Logical metabolite pattern

Anything you might ask an expert LC-MSMS reviewer to check.....

---

---

---

---


---

---

---

---

## Workflow detail

7. Field Mapping  = round, red thing = tomato

MS field name=LIS field name or HL7 segment

Instrument	HL7	LIS
"Sample ID"	= "PID-2-1"	= "Accession Number"
"Serial Number"	= "OBR-31-2"	= "Instrument ID"
"Compound"	= "OBR-31-5"	= "Test Mnemonic"

---

---

---

---

---

---

---

---

## Workflow detail

8. Test transmission from MS to LIS/MW – troubleshoot & repeat

- Create test patient orders in the TEST LIS
- Do an LC-MS run with the test patient acc#s
- Review the data with the MS software
- Transmit through the interface to the TEST LIS
- Look at communication trace(s) (if direct to LIS)
- Review in the LIS
- Find the problems and fix them

---

---

---

---

---

---

---

---

## Workflow Detail

9, 10, 11. Setup in LIVE, retest (document), go live!

- The myth of “copy to LIVE”
- Any interface build with manual entry (configuration, rules, test names) needs to be tested in LIVE
- Does an LIS TEST system that matches the LIS LIVE system exactly exist anywhere?
- After Go Live - verify some production patients in LIS & EMR LIVE environments

---

---

---

---

---

---

---

---

## Liquid Handler Interfacing (? Orders Download ?)

- Use automated barcode readers
- Is import (order download from LIS/MW) needed?
- Avoid LIS worklists/batches that require samples to be sorted (by a human)
- Can meta-data from other instruments (e.g., drug immunoassay numeric values) be downloaded to control dilutions/extraction protocols on the liquid handler?

---

---

---

---

---

---


---

---



## Liquid Handler to MS interface

- Export the plate map from the liquid handler (.txt, .csv, .xls)
  - Well ID (A1, H10)
  - Accession #
  - Plate barcode
  - Test name
  - ? Liquid handler ID/destination LC-MSMS ID?
- Export liquid handling errors (clot detection, short sample)?
- Options:
  - import directly to the MS (copy/paste to sample table)
  - OR - Import to the LIS/MW (do LIS/MW drivers exist for your liquid handler?), reformat, and then transfer to the MS



---

---

---

---

---

---

---

---



---

---

---

---

---

---

---

---

## Automating Data Review

- Defining "LC-MSMS Data Review"
- Automating data review & transfer with Rules Engine(s) in the:
  - LC-MSMS software
  - 3<sup>rd</sup> party software
  - MW
  - LIS

---

---

---

---


---

---

---

---

## 1st - Peak Review




---

---

---

---

---

---


---

---

## Using LC-MSMS metadata #1

### • MRM qualifier/quantifier ratio

Quantifier MRM (337/105) peak area = 2.8 e6  
 Qualifier MRM (337/188) peak area = 1.2 e6  
 Qual Peak Area ÷ Quant Peak Area Ratio = 0.43  
 Range = Expected +/- (20% to 30%)



### • Relative Retention Time

Retention Time Analyte ÷ Retention Time of Internal Standard  
 For example: lorazepam Rt 3.47 min ÷  
 lorazepam-d4 Rt 3.45 min = 1.01  
 Range = Expected +/- (0.02 min or 3%)

### • Expected = Average of MRM Ratios or RRTs found in calibrators from each run

---

---

---

---

---


---

---

---

## Using LC-MSMS metadata #2

- Use Internal Standard (I.S.) Peak Areas as an indicator of ion suppression, pipetting errors (no I.S.; twice the I.S.), injection errors (no I.S.)
- Expected I.S. Peak Areas = average of calibrator I.S. Peak Areas in the batch
- Acceptable range = derive from validation data and consensus documents (e.g. +/- 50%)



CAP CHM.18850  
 and CHM.18900

---

---

---

---

---

---

---

---

### Using LC-MSMS metadata #3

Additional parameters to consider:

- Abnormal peak shape: Peak Width, Peak Asymmetry
- Calibration errors:  $R^2$ , % Deviation (% Accuracy), slope, intercept
- Presence of Metabolites (e.g. EDDP present if Methadone present)
- Any quality parameter you ask an expert LC-MSMS reviewer to evaluate

---

---

---

---

---

---

---

### Setting targets & ranges

- Use “knowns” in each batch or across batches – i.e. calibrators – to set target means or thresholds for each batch
  - Expected RRt or Rt Rt 2.37 and RRt 1.01
  - Expected Qual MRM Ratio 43%
  - Expected I.S. Peak Area 2.34 e5
- Derive acceptable ranges from validation data and consensus standards (e.g. CLSI, WADA)
  - Rt +/- 0.2 min or +/- 3%; RRt +/- 0.02
  - Qual MRM Ratio +/- 25%
  - I.S. Peak Area +/- 50%

---

---

---

---

---

---

---

### Using rules engines – to do what?

1. Identify “quality failures” of LC-MSMS metadata
2. Prevent reporting of results with quality failures (or results outside of the AMR or for gender/age constraints) by:
  - a. **Don't export from the LC-MSMS (LC-MSMS vendor software rule)**
  - b. **OR don't transfer from MW to LIS (MW rule)**
  - c. **OR don't post results to the LIS patient record (LIS rule)**
  - d. **OR don't autoverify results in the LIS (LIS rule)**

---

---

---

---

---

---

---

## Rules engines – how to choose?

- Vendor software – minimal to extensive capability, optimal for quality failures, User friendly? User accessible?
- 3<sup>rd</sup> party generic software – highly flexible, IT resources?
- Middleware – potentially the most robust rules engine & more user accessible than vendor software or LIS but are MS data fields definable?
- LIS – usually least capability
- Commercial 3<sup>rd</sup> party MS data review automation – most sophisticated, but user accessible? \$\$\$

---

---

---

---

---

---

---

---

---

---

## 3<sup>rd</sup> party SW to ID quality failures




Fig. 4. Flow chart defining the quality control algorithm.

\* Dickerson JA, Schmeling M, Hoofnagle AN, Hoffman NG; Design and implementation of software for automated quality control and data analysis for a complex LC/MS/MS assay for urine opiates and metabolites, Clin Chim Acta , 415 (2013) 290–294.

---

---

---

---

---

---

---

---

---

---

## LC-MSMS Vendor SW to ID quality failures

Name	Sample	Concentration	Area	Failure Code	Used
100	100	100	100	0	
101	101	101	101	0	
102	102	102	102	0	
103	103	103	103	0	
104	104	104	104	0	
105	105	105	105	0	
106	106	106	106	0	
107	107	107	107	0	
108	108	108	108	0	
109	109	109	109	0	
110	110	110	110	0	
111	111	111	111	0	
112	112	112	112	0	
113	113	113	113	0	
114	114	114	114	0	
115	115	115	115	0	
116	116	116	116	0	
117	117	117	117	0	
118	118	118	118	0	
119	119	119	119	0	
120	120	120	120	0	
121	121	121	121	0	
122	122	122	122	0	
123	123	123	123	0	
124	124	124	124	0	
125	125	125	125	0	
126	126	126	126	0	
127	127	127	127	0	
128	128	128	128	0	
129	129	129	129	0	
130	130	130	130	0	
131	131	131	131	0	
132	132	132	132	0	
133	133	133	133	0	
134	134	134	134	0	
135	135	135	135	0	
136	136	136	136	0	
137	137	137	137	0	
138	138	138	138	0	
139	139	139	139	0	
140	140	140	140	0	
141	141	141	141	0	
142	142	142	142	0	
143	143	143	143	0	
144	144	144	144	0	
145	145	145	145	0	
146	146	146	146	0	
147	147	147	147	0	
148	148	148	148	0	
149	149	149	149	0	
150	150	150	150	0	
151	151	151	151	0	
152	152	152	152	0	
153	153	153	153	0	
154	154	154	154	0	
155	155	155	155	0	
156	156	156	156	0	
157	157	157	157	0	
158	158	158	158	0	
159	159	159	159	0	
160	160	160	160	0	
161	161	161	161	0	
162	162	162	162	0	
163	163	163	163	0	
164	164	164	164	0	
165	165	165	165	0	
166	166	166	166	0	
167	167	167	167	0	
168	168	168	168	0	
169	169	169	169	0	
170	170	170	170	0	
171	171	171	171	0	
172	172	172	172	0	
173	173	173	173	0	
174	174	174	174	0	
175	175	175	175	0	
176	176	176	176	0	
177	177	177	177	0	
178	178	178	178	0	
179	179	179	179	0	
180	180	180	180	0	
181	181	181	181	0	
182	182	182	182	0	
183	183	183	183	0	
184	184	184	184	0	
185	185	185	185	0	
186	186	186	186	0	
187	187	187	187	0	
188	188	188	188	0	
189	189	189	189	0	
190	190	190	190	0	
191	191	191	191	0	
192	192	192	192	0	
193	193	193	193	0	
194	194	194	194	0	
195	195	195	195	0	
196	196	196	196	0	
197	197	197	197	0	
198	198	198	198	0	
199	199	199	199	0	
200	200	200	200	0	

---

---

---

---

---

---

---

---

---

---

## MW rules to prevent Tx to LIS

### If → Then

#### Rules Engine Syntax


If {Test Result of} "VD2" >"5" AND {MRMR} =  
{"FAIL"} Then {Suppress that test}  
AND {Print Failure Log}

#### Description

"If the 25-OH-Vitamin D2 result is >LLOQ and there is an MRM Ratio Failure – then do not post the result in the LIS and print an error log"

## Automating Data Review with LIS rules

- Waters MassLynx HL7 Interface & SCC SOFT LIS
- Use MassLynx flagging functionality to identify MRM ratio, RRT, Internal Standard Peak area, R<sup>2</sup>, No Solution failures
- Send test & flag results to SOFT LIS across HL7 interface (customized data export)
- Write autoverification rules in SOFT LIS. Results with flagging failures inappropriate for concentration don't autoverify.




25-OH-D2

Low IS	1° Ratio Flag	R.T. Tolerance Flag	CD Flag	No Solution Flag
NO	YES	NO	NO	NO

Entire batch transmitted – problem samples must be suppressed – an individual sample **CANNOT** be removed or excluded from transmission.

## Commercial solutions for automating LC-MSMS data review


Indigo Biosystems ASCENT (cloud based)



## Workflow detail - #11to #14

Test Rules, Fix Rules & Retest (document)

- Write & follow a testing plan
- First test only **one** rule per result, then test **multiple** rules per result
- After edits – also retest what was working before (**regression** testing)
- Test batches of **genuine** **patient** samples that include “real life” rule failures



## Workflow detail - #11to #14

Test batches, scale up (document)

- What works for 1-10 samples may not work for 96 samples
- The best test batch for 2<sup>nd</sup> tier testing is a typical patient batch with production acc#s replaced by test acc#s

### Failures of interface testing -1

1. 5 MS error flags for autoverification rules
2. Tested one error per sample - YES Tested multiple errors per sample - NO
3. "No Peak (NP)" error flag rule – **autoverifies < LLOQ** as the result when no peak is present (no result sent by instrument but NP flag sent as "YES")
4. "Low Internal Standard (I.S.)" error flag rule – **result NOT autoverified** when I.S. is < threshold (I.S. flag sent as "YES")
5. When **both** I.S. and NP flags = YES – result <LLOQ should not autoverify but does because NP flag rule fires first (rule hierarchy). **A sample with no internal standard peak was incorrectly reported as <LLOQ.**

---

---

---

---

---

---

---

---

### Failures of interface testing - 2

- 10 to 20 samples/batch tested successfully in both TEST and LIVE environments
- Ordering test patients restricted to LIS personnel – no LIS resources to create 96, 192, 288 test orders
- Production starts – every 2-5 batches of 96 fails transmission to middleware (driver code fix needed)

---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

## Summary

- Productivity↑ can also = Error↑
- Write a testing plan – Test, Test, Test – Document, Document, Document
- No one knows your MS work flow better than you – why not be involved in customizing the automation of that workflow to get exactly what you want?
- You can do it!!

---

---

---

---

---

---

---